

Product Datasheet

Swine IL-1 alpha Recombinant Protein (orb1820389)

Description	The Swine IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SATYSFQSNM KYNFMRVINH QCILNDARNQ SIIRDPSGQY LMAAVLNNLD EAVKFDMAAY TSNDDSQLPV TLRISETRLF VSAQNEDEPV LLKELPETPK TIKDETSLLF FWEKHGNMDY FKSAHPKLF IATRQEKL VH MAPGLPSVTD FQILENQS (158)
Purity	98%
Source	Yeast
MW	18.1 kDa
Entrez	397094
Expiration Date	6 months from date of receipt.