

Product Datasheet

Zebrafish IL-4 Recombinant Protein (orb1820394)

Description	The Zebrafish IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SKLKIDEILL MEIIQSVNGI LNGKGEKMDL DQFIPDIYET GHYSKKTLCQ AGMALKGKIP SRSHLQRQLN AYAAYTGTGN CSVSTSGECT MKEFLEKTKA CCHYLYSAQR T (111)
Purity	98%
Source	Yeast
MW	12.4 kDa
Entrez	100188894
Expiration Date	6 months from date of receipt.