

Product Datasheet

Feline IL-17A Recombinant Protein (orb1820396)

Description	The Feline IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIAFPQNPGC PTTEDKNFPQ HVKVNINILN GNKSSRRPLD YYRRSTSPWS LHRNEDPERY PSVIWEAKCL HWGCVNTEGK EDHMMNSVPI QQEILVLRRE SRHCPHSFRL EKMLVTVGCT CVTPIVRHVV (130)
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	101095339
Expiration Date	6 months from date of receipt.