

Product Datasheet

Cynomolgus Monkey IL-17A Recombinant Protein (orb1820400)

Description	The Cynomolgus Monkey IL-17A yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIAIPRNSGC PNPSEDKNFPR TVMVNLNIHN RNTSTNPKRS SDYYNRSTSP WNLHRNEDPE RYPVSIWEAK CRHLGCVKAD GNVDYHMNSV PIQQEILVLR REPRHCPNSF RLEKILVSVG CTCVTPIVHH VA (132)
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	102119976
Expiration Date	6 months from date of receipt.