

Product Datasheet

Caprine IL-13 Recombinant Protein (orb1820412)

Description	The Caprine IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPSPVPSPSA LKELIEELVN ITQNQVPLCN GSMVWSLNLT SSMYCAALDS LISINCSVI QRTKRMLSAL CPHKPSAKHV SSEYIRDTKI EVAQFLKDLL RHRIVFRNE SYN (113)
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	102187477
Expiration Date	6 months from date of receipt.