

Product Datasheet

Caprine IL-4 Recombinant Protein (orb1820414)

Description	The Caprine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKCDITLLEEI IKMLNILTSQ KNSCMELPVA DVFAAPKNAT EKETFCRAGI ELRRIYRNHM CLNKFLGGLD RNLSSLASKT CSVNEAKTST STLRDLLERL KTIMKEKYSKC (111)
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	100860814
Expiration Date	6 months from date of receipt.