

Product Datasheet

Rabbit CCL3 Recombinant Protein (orb1820417)

Description	The Rabbit CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IPLGADTPTA CCFSYISRQI PYKFIADYFE TSSQCSKPGV IFLTKRGRQV CADISEAWVQ GYINDLELNS (70)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	100348776
Expiration Date	6 months from date of receipt.