

Product Datasheet

Swine LIF Recombinant Protein (orb1820424)

Description	The Swine LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLSITPVNA TCATRHPCHS NLMNQIKNQL AQLNSSANAL FILYYTAQGE PFPNNLDKLC GPNVTNFPPF HANGTEKARL VELYRIIAYL GASLGNITRD QRSLNPGAVN LHSKLNATAD SMRGLLSNVL CRLCNKYHVA HVDVAYGPDY SGKDVVFQKKK LGCQLLGKYK QVISVLARAF (180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	399503
Expiration Date	6 months from date of receipt.