

Product Datasheet

Rabbit LIF Recombinant Protein (orb1820426)

Description	The Rabbit LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPINPVNA TCNTHHPCPS NLMSQIRSQL AQLNGTANAL FILYYTAQGE PFPNNLDKLC GPNVTFDFPPF HANGTEKVRL VELYRIVAYL GTALGNITRD QKTLNPTAHS LHSKLNATAD TLRGLLSNVL CRLCSKYHVA HVDVAYGPD SGKDVVFQKKK LGCQLLGKYK QVMAVLAQAF (180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	100358914
Expiration Date	6 months from date of receipt.