

Product Datasheet

Canine LIF Recombinant Protein (orb1820428)

Description	The Canine Leukemia Inhibitory Factor yeast-derived recombinant protein is not tagged and is natural...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCHS NLMNQIRNQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFFPPF HANSTEKTRL VELYRIIAYL GASLGNITRD QKVLNPNALS LHSKLNATAD IMRGLLSNVL CRLCSKYHVA HVDVAYGPDT SGKDVVFQKKK LGCQLLGKYK QVIAVVAQAF (180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	403449
Expiration Date	6 months from date of receipt.