

## Product Datasheet

# Mouse IFN alpha Recombinant Protein (orb1820432)

<b>Description</b>	The Mouse IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
<b>Target</b>	IFN alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	CDLPQTHNLR NKRALLLVQ MRRLSPLSCL KDRKDFGFPO EKVDAQQIKK AQAIPVLSL TQQILNIFTS KDSSAAWNTT LLDSFCNDLH QQLNDLQGCL MQQVGVQEFP LTQEDALLAV RKYFHRITVY LREKKHSPCA WEVVRAEVWR ALSSSANVLG RLREEK (166)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.1 kDa
<b>Entrez</b>	<a href="#">15962</a>
<b>Expiration Date</b>	6 months from date of receipt.