

Product Datasheet

Mouse IFN alpha Recombinant Protein (orb1820433)

| | |
|-------------------------|---|
| Description | The Mouse IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free.... |
| Target | IFN alpha |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | CDLPQTHNLR NKRALLLVQ MRRLSPLSCL KDRKDFGFPO EKVDAQQIKK AQAIPVLSL TQQILNIFTS KDSSAAWNTT LLDSFCNDLH QQLNDLQGCL MQQVGVQEFP LTQEDALLAV RKYFHRITVY LREKKHSPCA WEVVRAEVWR ALSSSANVLG RLREEK (166) |
| Purity | 98% |
| Source | Yeast |
| MW | 19.1 kDa |
| Entrez | 15962 |
| Expiration Date | 6 months from date of receipt. |