

Product Datasheet

Mouse IFN alpha Recombinant Protein (orb1820434)

Description	The Mouse IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CDLPQTHNLR NKRALLLVQ MRRLSPLSCL KDRKDFGFPO EKVDAQQIKK AQAIPVLSL TQQILNIFTS KDSSAAWNTT LLDSFCNDLH QQLNDLQGCL MQQVGVQEFP LTQEDALLAV RKYFHRITVY LREKKHSPCA WEVVRAEVWR ALSSSANVLG RLREEK (166)
Purity	98%
Source	Yeast
MW	19.1 kDa
Entrez	15962
Expiration Date	6 months from date of receipt.