

Product Datasheet

Bovine LIF Recombinant Protein (orb1820435)

Description	The Bovine Leukemia Inhibitory Factor yeast-derived recombinant protein is not tagged and is natural...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCPS NLMNQIRNQL GQLNSSANSL FILYYTAQGE PFPNNLDKLC SPNVTDFFPPF HANATEKARL VELYRIIAYL GASLGNITRD QKVLNPYAHG LHSKLSTTAD VLRGLLSNVL CRLCSKYHVS HVDVTYGPDT SGKDVFQKKK LGCQLLGKYK QVIAVLAQAF (180)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	280840
Expiration Date	6 months from date of receipt.