

Product Datasheet

Canine IL-4 Recombinant Protein (orb1820439)

Description	The Canine IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HNFNITIKEI IKMLNILTAR NDSCMELTVK DVFTAPKNTS DKEIFCRAAT VLRQIYTHNC SNRYLRGLYR NLSSMANKTC SMNEIKKSTL KDFLERLKVI MQKKYYRH (108)
Purity	98%
Source	Yeast
MW	12.8 kDa
Entrez	403785
Expiration Date	6 months from date of receipt.