

Product Datasheet

Rabbit CD40 Ligand Recombinant Protein (orb1820443)

Description	The Rabbit CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	CD40 Ligand
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	<pre>MQKGDQDPQI AAHLISEASS KSSSVLQWAK KGYTMSNTL VTLENGKQLK VKRQGFYIY AQVTFCSNQE PSSQAPFIAS LCLKSSGGSE RILLRAANAR SSSKTCEQQS IHLGGVFELQ ADASVFNVT DASQVNHGTG FTSFGLLKL (149)</pre>
Purity	98%
Source	Yeast
MW	16.2 kDa
Entrez	100358388
Expiration Date	6 months from date of receipt.