

## Product Datasheet

# Rabbit CD40 Ligand Recombinant Protein (orb1820444)

<b>Description</b>	The Rabbit CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
<b>Target</b>	CD40 Ligand
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	<pre>MQKGDQDPQI AAHLISEASS KSSSVLQWAK KGYTMSNTL VTLENGKQLK VKRQGFYIY AQVTFCSNQE PSSQAPFIAS LCLKSSGGSE RILLRAANAR SSSKTCEQQS IHLGGVFELQ ADASVFNVT DASQVNHGTG FTSFGLLKL (149)</pre>
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.2 kDa
<b>Entrez</b>	<a href="#">100358388</a>
<b>Expiration Date</b>	6 months from date of receipt.