

## Product Datasheet

# Human/Monkey CD40 Ligand (CD40LG) Recombinant Protein (orb1820447)

<b>Description</b>	The Human & Monkey CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endo...
<b>Target</b>	CD40LG
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MQKGDQNPQI AAHWISEASS KTTSVLQWAE KGYTMSNNL VTLENGKQLT VKRQGLYIY AQVTFCSNRE ASSQAPFIAS LCLKSPGRFE RILLRAANTH SSAKPCGQQS IHLGGVFELQ PGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.2 kDa
<b>Entrez</b>	<a href="#">959</a>
<b>Expiration Date</b>	6 months from date of receipt.