

## Product Datasheet

# Feline CD40 Ligand (CD40LG) Recombinant Protein (orb1820448)

|                         |                                                                                                                                                                                              |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | The Feline CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...                                                                                      |
| <b>Target</b>           | CD40 Ligand                                                                                                                                                                                  |
| <b>Form/Appearance</b>  | Lyophilized                                                                                                                                                                                  |
| <b>Storage</b>          | -20°C                                                                                                                                                                                        |
| <b>Note</b>             | For research use only                                                                                                                                                                        |
| <b>Protein Sequence</b> | MQKGDQDPRV AAHVISEASS<br>STASVLQWAP KGYTTISNL<br>VTLENGKQLA VKRQGLYYIY<br>AQVTFCSNRE ASSQAPFIAS<br>LCLHSPSGSE RVLLRAANAR<br>SSSKPCGQQS IHLGGVFELH<br>PGASVFNVT DPSQVSHGTG<br>FTSFGLLKL (149) |
| <b>Purity</b>           | 98%                                                                                                                                                                                          |
| <b>Source</b>           | Yeast                                                                                                                                                                                        |
| <b>MW</b>               | 15.9 kDa                                                                                                                                                                                     |
| <b>Entrez</b>           | <a href="#">493850</a>                                                                                                                                                                       |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                                                                                                               |