

## Product Datasheet

# Equine CD40 Ligand (CD40LG) Recombinant Protein (orb1820449)

|                         |   |
|-------------------------|---|
| <b>Description</b>      | The Equine CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...   |
| <b>Target</b>           | CD40LG  |
| <b>Form/Appearance</b>  | Lyophilized   |
| <b>Storage</b>          | -20°C   |
| <b>Note</b>             | For research use only   |
| <b>Protein Sequence</b> | MQKGDQEPQI AAHVISEASS<br>KTASVLQWAQ KGYTTISNNL<br>VTLENGKQLA VKRQGLYYIY<br>AQVTFCSNRE ASGQAPFIAS<br>LCLRSVSGSE RILLRAANTH<br>SSSKPCGQQS IHLGGVFELQ<br>PGASVFNVT DPSQVSHGTG<br>FTSFGLLKL (149) |
| <b>Purity</b>           | 98%   |
| <b>Source</b>           | Yeast   |
| <b>MW</b>               | 16.0 kDa  |
| <b>Entrez</b>           | <a href="#">100054624</a>   |
| <b>Expiration Date</b>  | 6 months from date of receipt.  |