

Product Datasheet

Chicken CX3CL1 Recombinant Protein (orb1820450)

Description	The Chicken CX3CL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	CX3CL1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QPRAPLKCSK WCISFHRAID QRQIKSYRET EPQCTKKAI FTTKRNREIC ANPYEPWVEK IVKKLD (66)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	415651
Expiration Date	6 months from date of receipt.