

Product Datasheet

Human CCL5 Recombinant Protein (orb1820453)

Description	The Human CCL5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	CCL5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPYSSDTTPC CFAYIARPLP RAHIKEYFYT SGKCSNPAVV FVTRKNRQVC ANPEKKWVRE YINSLEMS (68)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	6352
Expiration Date	6 months from date of receipt.