

Product Datasheet

Canine CCL3 Recombinant Protein (orb1820458)

Description	The Canine CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QVFSPFGADT PIACCFYSYVS KQIPRKFIVD CFETSSQCCK PGIIFETRKG RQACANPSEA WVQEYVADLK LKA (73)
Purity	98%
Source	Yeast
MW	8.1 kDa
Entrez	448787
Expiration Date	6 months from date of receipt.