

## Product Datasheet

### Dolphin IL-10 Recombinant Protein (orb1820460)

<b>Description</b>	The Dolphin IL-10 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
<b>Target</b>	IL-10
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GQHQGTQAKD SCIHFPDSL HMLRELRAAF SSVKTFQMN DQLDNSLLSQ SLLEDFKGYL GCQALSEMIQ FYLEEVMPKA EDHGPNKEH VNLSGEKLT LRLQLRRCHR FLPCENKSKA VEQVKSVMFNK LQKEGVYKAM SEFDIFINYI EAYMITKMKN (160)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.6 kDa
<b>Entrez</b>	<b>101330184</b>
<b>Expiration Date</b>	6 months from date of receipt.