

Product Datasheet

Dolphin IL-10 Recombinant Protein (orb1820461)

Description	The Dolphin IL-10 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-10
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GQHQGTQAKD SCIHFPDSLP HMLRELRAAF SSVKTFQMN DQLDNSLLSQ SLLEDFKGYL GCQALSEMIQ FYLEEVMPKA EDHGPNIEH VNLSLGEKLT LRLQLRRCHR FLPCENKSKA VEQVKSVMFNK LQKEGVYKAM SEFDIFINYI EAYMITKMKN (160)
Purity	98%
Source	Yeast
MW	18.6 kDa
Entrez	101330184
Expiration Date	6 months from date of receipt.