

Product Datasheet

Bovine CCL3 Recombinant Protein (orb1820474)

Description	The Bovine CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APFGADTPTA CCFSYVARQL SRKIVADYFE TSSQCSKPGV IFQTKKGRQV CANPTEDWVQ EYITDLELNA (70)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	282170
Expiration Date	6 months from date of receipt.