

## Product Datasheet

# Chicken CD40 Ligand (CD40LG) Recombinant Protein (orb1820480)

<b>Description</b>	The Chicken CD40 Ligand Ligand yeast-derived recombinant protein is not tagged and is naturally endo...
<b>Target</b>	CD40LG
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MHRGHEHPHL KSRNETSVAE EKRQPIATHL AGVKSNTTVR VLKWMTTSYA PTSSLISYHE GKLKVEKAGL YYIYSQVSFC TKAAASAPFT LYIYLYLPME EDRLLMKGLD THSTSTALCE LQSIREGGVF ELRQGDMVFV NVDSTAVNV NPGNTYFGMF KL (162)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.2 kDa
<b>Entrez</b>	<a href="#">395485</a>
<b>Expiration Date</b>	6 months from date of receipt.