

Product Datasheet

Chicken CD40 Ligand (CD40LG) Recombinant Protein (orb1820482)

Description	The Chicken CD40 Ligand Ligand yeast-derived recombinant protein is not tagged and is naturally endo...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MHRGHEHPHL KSRNETSVAE EKRQPIATHL AGVKSNTTVR VLKWMTTSYA PTSSLISYHE GKLKVEKAGL YYIYSQVSFC TKAAASAPFT LYIYLYLPME EDRLLMKGLD THSTSTALCE LQSIREGGVF ELRQGDMVFV NVDSTAVNV NPGNTYFGMF KL (162)
Purity	98%
Source	Yeast
MW	18.2 kDa
Entrez	395485
Expiration Date	6 months from date of receipt.