

Product Datasheet

Canine CD40 Ligand (CD40LG) Recombinant Protein (orb1820485)

Description	The Canine CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MQKGDQDPRI AAHWISEASS NPASVLRWAP KGYTTISSNL VSLENGKQLA VKRQGLYYVY AQVTFCSNRA ASSQAPFVAS LCLHSPSGTE RVLLRAASSR GSSKPCGQQS IHLGGVFELH PGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	15.8 kDa
Entrez	403468
Expiration Date	6 months from date of receipt.