

Product Datasheet

Swine IL-2 Recombinant Protein (orb1820489)

Description	The Swine IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTKNT KKQLEPLLLD LQLLLKEVKN YENADLSRML TFKFYMPKQA TELKHLQCLV EELKALEGVL NLGQSKNSDS ANIKESMNNI NVTVLELKGS ETSFKCEYDD ETVTAVEFLN KWITFCQSIY STL T (134)
Purity	98%
Source	Yeast
MW	15.2 kDa
Entrez	396868
Expiration Date	6 months from date of receipt.