

Product Datasheet

Swine CD40 Ligand (CD40LG) Recombinant Protein (orb1820495)

Description	The Swine CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MHKGDQDPQI AAHWISEASS KTASVLQWAP KGYTTLSTNL VTLENGRQLA VKRQGIYYIY AQVTFCSNRD AAGQAPFIAS LCLRSPSGSE RILLRAANTH SSSKPCGQQS IHLGGVFELQ PGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.0 kDa
Entrez	397231
Expiration Date	6 months from date of receipt.