

Product Datasheet

Canine Erythropoietin Recombinant Protein (orb1820504)

Description	The Canine Erythropoietin yeast-derived recombinant protein is not tagged and is naturally endotoxin...
Target	EPO
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APRLICDSR VLERYILEAR EAENVTMGCA QGCSFSENIT VPDTKVN FYT WKRMDVGQQA LEVWQGLALL SEAILRGQAL LANASQPSET PQLHVDKAVS SLRSLTLLR ALGAQKEAMS LPEEASPAPL RTFTVDTLCK LFRIYSNFLR GKLTLYTGEA CRRGDR (166)
Purity	98%
Source	Yeast
MW	18.4 kDa
Entrez	404002
Expiration Date	6 months from date of receipt.