

Product Datasheet

Human IL-17A Recombinant Protein (orb1820519)

Description	The Human IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GITIPRNPGC PNSEDKNFPR TVMVNLNIHN RNTNTNPKRS SDYYNRSTSP WNLHRNEDPE RYPSVIWEAK CRHLGCINAD GNVDYHMNSV PIQQEILVLR REPPHCPNSF RLEKILVSVG CTCVTPIVHH VA (132)
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	3605
Expiration Date	6 months from date of receipt.