

Product Datasheet

Canine IL-17A Recombinant Protein (orb1820522)

Description	The Canine IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPQNP GCRNT EDKNFPQHVK VNLNILNRNT NSRRPSDYNN RSTSPWNLHR NEDPERYPSV IWEAKCRHLG CVNNEGNINY HMNSVPIQQE ILVLRRESQH CPHSFRLEKM LVAVGCTCVT PIVRHVA (127)
Purity	98%
Source	Yeast
MW	14.8 kDa
Entrez	481837
Expiration Date	6 months from date of receipt.