

Product Datasheet

Canine GM-CSF Recombinant Protein (orb1820530)

Description	The Canine GM-CSF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	GM-CSF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTRSPTLVT RPSQHVDIAQ EALSLNNSN DVTAVMNKAV KVVSEVFDPE GPTCLETRLQ LYKEGLQGSL TSLKNPLTMM ANHYSQHCPP TPESPCATQN INFKSFKENL KDFLFNIPFD CWKPVKK (127)
Purity	98%
Source	Yeast
MW	14.2 kDa
Entrez	403923
Expiration Date	6 months from date of receipt.