

Product Datasheet

Mouse IL-1 alpha Recombinant Protein (orb1820532)

Description	The Mouse IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SAPYTYQSDL RYKLMKLVQR KFVMNDSLNQ TIYQDVVKHY LSTTWLNDLQ QEVKFDYAY SSGGDDSKYP VTLKISDSQL FVSAQGEDQP VLLKELPETP KLITGSETDL IFFWKSINSK NYFTSAAYPE LFIATKEQSR VHLARGLPSM TDFQIS (156)
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	16175
Expiration Date	6 months from date of receipt.