

## Product Datasheet

# Mouse IL-1 alpha Recombinant Protein (orb1820533)

<b>Description</b>	The Mouse IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IL-1 alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SAPYTYQSDL RYKLMKLVQR KFVMNDSLQ TIYQDVKHY LSTTWLNDLQ QEVKFDYAY SSGGDDSKYP VTLKISDSQL FVSAQGEDQP VLLKELPETP KLITGSETDL IFFWKSINSK NYFTSAAYPE LFIATKEQSR VHLARGLPSM TDFQIS (156)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	18.0 kDa
<b>Entrez</b>	<a href="#">16175</a>
<b>Expiration Date</b>	6 months from date of receipt.