

Product Datasheet

Ovine TNF alpha Recombinant Protein (orb1820541)

Description	The Ovine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LRSSSQASNN KPVAVHVANI SAPGQLRWGD SYANALMANG VELKDNQLVV PTDGLYLIYS QVLFRGHGCP STPLFLTHTI SRIAVSYQTK VNILSAIKSP CHRETLEGAE AKPWYEPIYQ GGVFQLEKGD RLSAEINLPE YLDYAESGQV YFGIIAL (157)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	443540
Expiration Date	6 months from date of receipt.