

Product Datasheet

Mouse TWEAK Recombinant Protein (orb1820551)

| | |
|-------------------------|--|
| Description | The Mouse TWEAK yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | TWEAK |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SAPKGRKARP RRAIAAHYEV HPRPGQDGAQ AGVDGTVSGW EETKINSSSP LRYDRQIGEF TVIRAGLYYL YCQVHFDEGK AVYLKLDLLV NGVLALRCLE EFSATAASSP GPQLRLCQVS GLLPLRPGSS LRIRTLPSWAH LKAAPFLTYF GLFQVH (156) |
| Purity | 98% |
| Source | Yeast |
| MW | 17.2 kDa |
| Entrez | 21944 |
| Expiration Date | 6 months from date of receipt. |