

Product Datasheet

Swine IL-12 p40 Recombinant Protein (orb1820552)

Description	The Swine IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally endotoxin...
Target	IL-12 p40
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IWELEKNVYV VELDWYPNAP GEMVVLTCNT PEEDGITWTS DQSSEVLGTG KTLTIHVKEF GDAGQYTCRK GGAVLSQSL LLHKKEDGIW STDILKDQKE PKNKSFLKCE AKNYSGRFTC WWLTAISTDL KFSVKSSRGS TDPRGVTCGT ATLSIDLGEY KKYRVECQEG SACPAAEESL PIEVVLEAVH KLKYENYTSS FFIRDIKPD PPKNLQLNPL KNSRHVEISW EYPDTWSTPH SYFSLMFGVQ VQGKMKREKK DKLFTDQISA KVTCHKDANI RVQARDRYYS SSWSEWASVS CN (302)
Purity	98%
Source	Yeast
MW	34.3 kDa
Entrez	397076
Expiration Date	6 months from date of receipt.