

Product Datasheet

Mouse IL-18 Recombinant Protein (orb1820553)

Description	The Mouse IL-18 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-18
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NFGRLHCTTA VIRNINDQVL FVDKRQPVFE DMTDIDQSAS EPQTRLIYM YKDSEVRGLA VTLSVKDSKM STLSCKNKII SFEEMDPEN IDDIQSDLIF FQKRVPGHMK MEFESSLYEG HFLACQKEDD AFKLILKKKD ENGDKSVMFT LTNLHQS (157)
Purity	98%
Source	Yeast
MW	18.1 kDa
Entrez	16173
Expiration Date	6 months from date of receipt.