

Product Datasheet

Human IL-1 beta Recombinant Protein (orb1820554)

Description	The Human IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IL-1 beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APVRSLNCTL RDSQQKSLVM SGPYELKALH LQGQDMEQQV VFSMSFVQGE ESNDKIPVAL GLKEKNLYLS CVLKDDKPTL QLESVDPKNY PKKKMEKRFV FNKIEINNKL EFESAQFPNW YISTSQAENM PVFLGGTKGG QDITDFTMQF VSS (153)
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	3553
Expiration Date	6 months from date of receipt.