

Product Datasheet

Mouse IL-34 Recombinant Protein (orb1820557)

Description	The Mouse IL-34 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-34
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NENLEIWTLT QDKECDLTGY LRGKLQYKNR LQYMKHYFPI NYRIAVPYEG VLRVANITRL QKAHVSEREL RYLWVLVSLN ATESVMDVLL EGHPSWKYLQ EVQTLLENVQ RSLMDVEIGP HVEAVLSLLS TPGLSLKLVR PKALLDNCFR VMELLYCSCC KQSPILKWQD CELPRLPHPS PGSLMQCTAT NVYPLSRQTP TSLPGSPSSS HGSLP (215)
Purity	98%
Source	Yeast
MW	24.6 kDa
Entrez	76527
Expiration Date	6 months from date of receipt.