

Product Datasheet

Guinea Pig IL-15 Recombinant Protein (orb1820560)

Description	The Guinea Pig IL-15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-15
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SWHDVIYDLK RIENLIQSIH IDATLYTESD VHPNCKITAM KCFLELRVI LHESRNEDIH ETITNLILA NSSLNSNGNV TESGCKECEE LEEKSIAEFL QSFVHIVQMF INTS (114)
Purity	98%
Source	Yeast
MW	13.1 kDa
Entrez	100379533
Expiration Date	6 months from date of receipt.