

Product Datasheet

Human IL-6 Recombinant Protein (orb1820565)

Description	The Human IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VPPGEDSKDV AAPHRQPLTS SERIDKQIRY ILDGISALRK ETCNKS NMCE SSKEALAENN LNLPKMAEKD GCFQSGFN EE TCLVKIITGL LEFEVYLEYL QNRFESSEEQ ARAVQMSTKV LIQFLQKKAK NLDAITTPDP TTNASLLTKL QAQNQWLQDM TTHLILRSFK EFLQSSLRAL RQM (183)
Purity	98%
Source	Yeast
MW	20.8 kDa
Entrez	3569
Expiration Date	6 months from date of receipt.