

Product Datasheet

Human CCL4 Recombinant Protein (orb1820569)

| | |
|-------------------------|---|
| Description | The Human CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ... |
| Target | CCL4 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | APMGSDPPTA CCFSYTARKL PRNFVVDYYE TSSLCSQPAV VFQTKRSKQV CADPSETWVQ EYVYDLELN (69) |
| Purity | 98% |
| Source | Yeast |
| MW | 7.8 kDa |
| Entrez | 6351 |
| Expiration Date | 6 months from date of receipt. |