

Product Datasheet

Mouse CCL4 Recombinant Protein (orb1820574)

Description	The Mouse CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMGSDPPTS CCFSYTSRQL HRSFVMDYYE TSSLCSKPAV VFLTKRGRQI CANPSEPWVT EYMSDLELN (69)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	20303
Expiration Date	6 months from date of receipt.