

Product Datasheet

Rat CCL4 Recombinant Protein (orb1820576)

Description	The Rat CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The pu...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APIGSDPPTS CCFSYTSRKI HRNFVMDYYE TSSLCSQPAV VFLTKKGRQI CADPSEPWWN EYVNDLELN (69)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	116637
Expiration Date	6 months from date of receipt.