

Product Datasheet

Canine IL-13 Recombinant Protein (orb1820602)

Description	The Canine IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPSPVTPSPT LKELIEELVN ITQNQASLCN GSMVWSVNLT AGMYCAALES LINVSDCSAI QRTQRMLKAL CSQKPAAGIS SERSRDTKIE VIQLVKLLT YVRGVYRHGN FR (112) [Gln78/79Deletion]
Purity	98%
Source	Yeast
MW	12.3 kDa
Entrez	442990
Expiration Date	6 months from date of receipt.