

Product Datasheet

Bovine IL-17F Recombinant Protein (orb1820610)

Description	The Bovine IL-17F yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-17F
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QKPPRSGDSG LCPTLEDHIV RVDIRIRRN QGVSF SRNLQ NRSISPWDYN ITRDANRFPS EIAEAQCRHM GCINAEGQED RTLNSVPIQQ EFLVLR RKQN GCPGLFQLEK VLVTVGCTCV TPMVRHLHDP MADQVADQLS (140) [Phe34Ser]
Purity	98%
Source	Yeast
MW	15.9 kDa
Entrez	506030
Expiration Date	6 months from date of receipt.